

Measuring Instructions

LegAssist™ - LCS (Lobule Compression System)

STEP 3 - Measure Lengths (continued)

5. **Lateral Lengths (A2):** Measure lower leg from your mark on lateral aspect of Ø Point to your mark on lateral aspect of **Bottom of Garment** and record length on **Line A2**.
6. **Medial Length (B2):** Measure lower leg from your mark on the **medial** aspect of Ø Point to your mark on **medial** aspect of **Bottom of Garment** and record length on **Line B2**.
7. **Posterior Length (C2):** Measure lower leg from your mark on the **posterior** aspect of Ø Point to your mark on **posterior** aspect of **Bottom of Garment** and record length on **Line C2**.
8. **Anterior Lengths (D2):** Measure lower leg from your mark on **anterior** aspect of Ø Point to your mark on **anterior** aspect of **Bottom of Garment** and record length on **Line D2**.

STEP 4 - Measure Circumferences

1. **Upper:** Encircle leg with tape measure at each mark located on the **lateral** aspect between the Ø Point and **Gluteal Fold**, and record each circumference on the corresponding line, as applicable.
2. **Lower:** Encircle leg with tape measure at each mark located on the **lateral** aspect between the Ø Point and **Bottom of Garment**, and record each circumference on the corresponding line, as applicable.

STEP 5 - Hip attachment (Optional)

1. Measure from **lateral** aspect of **Gluteal Fold** to **lateral** aspect of hip at the desired location for the bottom of the belt and record length on **line F**.
2. Measure circumference of torso at the desired location for the belt and record circumference on **line G**.