



Call: (571) 232-3179

Fax: (703) 852-4371

Email: nancy@bodyworks-compression.com



Custom Measurement Form for Stump Shrinkers

Account Information

Account Number _____ Date _____

Account Name _____ Contact _____

Phone _____ Fax _____

Patient ID _____

Prescribing Physician _____

Circular Knit Compression Garment Information (without seam)

Compression: 20-30 mmHg Model: Dynamic (Varin Soft-in)
 30-40 mmHg Dynamic Silver (Varin Soft-in)


Style: C-F below knee F-G above knee H-D Pirogoff
 Silicone Border Hip Attachment

K-T Body Part (worn with F-G) Slip-on Hook and Loop Closure
 20-30 mmHg 30-40 mmHg

Extremity: Right Left Quantity: _____ Pieces

Comments: _____

Circumferences



Lengths

- S-T
- S-H
- S-K }
- S-G }
- S-F
- S-E
- S-D
- S-C
- S-B1
- S-B

Measure circumferences beginning 5 cm from the distal end of the stump.

S = Distal end of the stump. Measure all lengths from 'S'.