

P.O. #	Account #	Contact	Date
--------	-----------	---------	------

Phone	Patient Name	Age	Height	Weight
-------	--------------	-----	--------	--------

Company Name

Dropship Name & Address

PRODUCT INFORMATION

COMPREFLEX REDUCE TC	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
COMPREFLEX REDUCE KNEE	Size		Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

A PDF of this order form can be found online at bit.ly/SIGmccOrderforms

MEASURING INSTRUCTIONS


COMPREFLEX REDUCE THIGH COMPONENT

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the lateral aspect **Top of Patella** and continue to the 23cm dot. Record in the corresponding lines D and F (Left/Right). Measure length from the lateral aspect **Top of Patella** to lateral aspect **Gluteal Fold**, and record in box G1.

COMPREFLEX REDUCE KNEE

Measure circumferences where indicated by the black dots, beginning with the dot 5cm below the lateral aspect **Fibular Head**. Continue with the dot 5cm above the lateral aspect **Top of Patella**. Record in the corresponding lines C1 and D (Left/Right).

CIRCUMFERENCE



COMPREFLEX REDUCE TC SIZING CHART

	RANGE 1	RANGE 3
F	<80cm	<125cm
D	<65cm	<110cm

BLACK	LEFT		
	REGULAR	2401-TCR-L	2403-TCR-L
TALL	2401-TCT-L	2403-TCT-L	
RIGHT	REGULAR	2401-TCR-R	2403-TCR-R
	TALL	2401-TCT-R	2403-TCT-R
BEIGE	LEFT		
	REGULAR	2411-TCR-L	2413-TCR-L
TALL	2411-TCT-L	2413-TCT-L	
RIGHT	REGULAR	2411-TCR-R	2413-TCR-R
	TALL	2411-TCT-R	2413-TCT-R

COMPREFLEX REDUCE KNEE SIZING CHART

	RANGE 1	RANGE 3
D	<75cm	<115cm
C1	<65cm	<90cm
BLACK	2401-KP	2403-KP
BEIGE	2411-KP	2413-KP